

```
#####
## ChIP-seq Analysis Workflow ##
#####
## Mireia Ramos-Rodríguez    ##
## 14th of March 2019        ##
#####
```

```
#####
# Set up your working environment and necessary files
#####
# Open the virtual machine, log in to the VPN and connect to `ironwomen`.
# There, you will need to create the following directory structure.
```

```
#-----
# 1) Generate directory structure:
#-----
# yourpersonalfolder
# |-- ChIP-seq
#     |-- fastq
#     |-- BAM
#     |-- peaks
```

```
cd yourpersonalfolder
mkdir ChIP-seq
cd ChIP-seq
mkdir fastq
mkdir BAM
mkdir peaks
```

```
#-----
# 2) Copy sample file: NL1_h3k27ac.fastq.gz
#-----
cp ~/workshop_ChIPseq/ChIP-seq_sample/fastq/NL1_h3k27ac_sample.fastq.gz
~/yourpersonalfolder/ChIP-seq/fastq/.
```

```
#####
## Quality control with FastQC
#####
# First, we are going to check the quality of the sequenced reads using FastQC.
# The output is an html file with a summary of the tests performed and if your
# reads are passing or not the quality control thresholds.
```

```
#-----
# 1) Run FastQC on the sample file
#-----
cd ~/yourpersonalfolder/ChIP-seq # Make sure you are here
mkdir FastQC # Make directory for output
fastqc fastq/NL1_h3k27ac_sample.fastq.gz -o FastQC # 15 seconds
```

```
#-----
# 2) To see the results, copy the output html to the virtual machine
#-----
# Open a terminal in the virtual machine
scp msuser@ironwomen:/home/labs/mslab/msuser/yourpersonalfolder/ChIP-
seq/FastQC/NL1_h3k27ac_sample_fastqc.html .
# Now open the hmtl from your virtual machine and explore the results
```

```
#####
## Alignment with Bowtie2
#####
#-----
# 1) Align the sample file to the reference genome
#-----
bowtie2 -t -x ~/workshop_ChIPseq/reference_genome/hg19 -U fastq/NL1_h3k27ac_sample.fastq.gz
-S BAM/NL1_h3K27ac_sample.sam
```

```
#####
## Post-processing with Samtools
#####
#-----
# 1) Convert to BAM and sort
```

```

#-----
samtools view BAM/NL1_h3K27ac_sample.sam -b | samtools sort - -o
BAM/NL1_h3K27ac_sample.sorted.bam

#-----
# 2) Remove duplicates
#-----
markdup BAM/NL1_h3K27ac_sample.sorted.bam BAM/NL1_h3K27ac_sample.rmdup.bam -r -s
#-----
# 3) Index bam files
#-----
samtools index BAM/NL1_h3K27ac_sample.sorted.bam
samtools index BAM/NL1_h3K27ac_sample.rmdup.bam

#####
## Peak calling with MACS2
#####
activate-macs-git-2017.5.15 # Activate MACS2

mkdir tmp/ # Create folder for temporary files

#-----
# Run peak calling with MACS2
#-----
macs2 callpeak -f BAM -t BAM/NL1_h3K27ac_sample.rmdup.bam -c ~/workshop_ChIPseq/ChIP-
seq/BAM/NL1_input.sorted.bam -g hs -n peaks/NL1_h3K27ac --tempdir tmp/ --broad --nomodel

```